




UNICONE
PLUS




TRY NEW POSSIBILITIES.



Our product portfolio includes a wide range of rotary instruments and manual dental instruments.

We cover common and special needs of the dentists to make root canal treatment more simple.



FLEXIBILITY



CYCLIC FATIGUE RESISTANCE

NEW POSSIBILITIES THANKS TO NEW HEAT-TREATMENT OF NiTi.




reciprocal movement

- Complete system for reciprocal movement.
- Additional sizes.
- GP – new reciprocal file for glide path.



- New technology of heat treatment.
- Respect of the root canal anatomy.
- Safe to work with.




- Save time to reach working length.
- Cutting efficiency.

Unicone PLUS

COMPLETE RECIPROCAL SYSTEM FOR ROOT CANAL TREATMENT

Unicone GP

GLIDE PATH



 reciprocal
movement

taper ratio 4%
ISO 015




L = 25 144 50 ... 2815

Unicone PLUS

packing 6 pcs

SHAPING

FINISHING



 reciprocal
movement

taper ratio	6%	6%	6%	6%	4%	4%
ISO	020	025	035	040	035	040







L = 25 144 50 ... 2845 2855 2865 2875 2825 2835

Set III Unicone PLUS

MEDIN Endo Set

L = 25 397 144 45 1335

Sets contain 4 instruments




W-File; 010 L25	1x	
Unicone GP; 4/015	1x	
Unicone PLUS; 6/025	1x	
Unicone PLUS; 4/040	1x	

Set I Unicone PLUS

Unicone PLUS Set of 3

L = 25 397 144 45 1315

Sets contain 3 instruments




Unicone PLUS; 6/020	1x	
Unicone PLUS; 6/025	1x	
Unicone PLUS; 6/040	1x	

Set II Unicone PLUS

Unicone PLUS Set of 3

L = 25 397 144 45 1325

Sets contain 3 instruments

Unicone PLUS; 6/035	1x	
Unicone PLUS; 4/035	1x	
Unicone PLUS; 4/040	1x	

MEDIN, a.s.
Vlachovická 619
592 31 Nové Město na Moravě
Czech Republic

CZECH REPUBLIC
tel.: 566 684 248
fax: 566 684 384
e-mail: prodej@medin.cz

EXPORT
tel.: +420 566 684 336
fax: +420 566 684 385
e-mail: export@medin.cz

The logo for UNICONE PLUS features a stylized white cross with a small triangle at the top-left corner containing a plus sign. To the right of this symbol, the word "UNICONE" is written in a large, bold, white sans-serif font, and the word "PLUS" is written in a smaller, bold, white sans-serif font directly below it.

Distributor



www.medin.cz
www.medin.eu

